



# MiniMed Paradigm® Veo™ System

*A new era in diabetes management*



## The Pump that has it ALL !

- **EASY TO USE** – Very simple menu navigation and clear at-a-glance onscreen instructions
- **FLEXIBLE & PRECISE INSULIN DELIVERY** – puts you in control with low basal rate of 0.025 u/hr, programmable temporary basal rates and basal programmes to suit what you do everyday
- **HELPFUL** – with reminders and alerts to help even the busiest patients
- **SMART & SAFE** – Bolus Wizard® (insulin dose calculator) can help simplify bolus dosing & reduce errors

**DISCREET & FUN!** – small size with a large range of accessories, 5 different colours to choose from and a range of SKINS pump covers to be as individual as you are.

## Continuous Glucose Monitoring (CGM)

This is the real difference with a Medtronic Insulin Pump, and whether you choose to use CGM now or sometime in the future, the great news is that it you will already have it – built into your Medtronic Pump. To use whenever you like, as often as you like. CGM can allow you to take early action on rising or falling glucose levels\*, see how where and how quickly your levels are moving\*\*, and can alert you to take action when you reach levels that you want to avoid and minimize the risk of hyper- or hypoglycaemia<sup>1/2 \*\*\*</sup>

The Paradigm Veo System is the **First & Only** insulin pump available that is CGM-ready insulin pump, and features a totally unique automatic insulin shut-off mechanism !

## *Is Hypoglycaemia a Challenge?*

Many patients may still suffer from Hypoglycaemia Unawareness, Nocturnal Hypoglycaemia or a real fear of Hypoglycaemia. The Low Glucose Suspend (LGS) feature is designed to automatically detect hypoglycaemia and suspend insulin infusion prior to the onset of severe hypoglycaemia episodes. This is the exciting First component of the Closed Loop, and can really provide increased Peace of Mind for patients, parents and families.

## Medtronic Diabetes

- .....Medtronic supports Healthcare Professionals and their Patients
- .....Medtronic has a sound, secure and long History
- .....Medtronic provides Worldwide support
- .....Medtronic Invests in the Future and is Closing the Loop

Please visit [www.medtronic-diabetes.co.uk](http://www.medtronic-diabetes.co.uk)  
or contact us on 01923 205167 to see how we can help

\*Predictive alerts, \*\* Trend arrows, \*\*\*Low and High Threshold Alerts  
1.Garg S, Zisser H, Schwartz S. et. al., Improvements in glycaemic excursions with a transcutaneous, real-time continuous glucose sensor. *Diabetes Care*. 2006. 29:44-50. 2. Bode B, Gross K, Rikalo N, Schwartz S, Wahl T, Page C, Gross T, Mastrototaro J. Alarms Based on Real-Time Sensor Glucose Values Alert Patients to Hypo- and Hyperglycemia: The Guardian Continuous Monitoring System. *Diabetes Technol Ther*. 2004;6(2):105-113.

